

# In-cell UV cross-linking of photo-reactive ligands to RNA

---

Date: 7/20/2024

Author: Seth D. Veenbaas

Aim: Prepare ligand crosslinked total RNA for JuMP RT, and sequencing.

Logic: Treating cells with photo-reactive ligand and 365 nm UV cross-links ligand to RNA at direct binding sites, which can be subsequently mapped through JuMP RT and sequencing.

**NOTE: All experiments should start fresh by streaking B. sub on plates.**

- B. sub does not grow well from glycerol stocks or passaged cell pellets.
- B. sub does not store well at 4 °C because it enters a spore form.

## Day 1: Plate cells

Materials:

- B. subtilis str. 168 glycerol stock
- LB Agar plate

Protocol:

1. Streak a plate with B. subtilis str. 168 to form single colonies.
2. Incubate plate at 30 °C overnight.

Notes:

## Day 2: Overnight culture

Materials:

- LB media (5 mL)
- Culture tube

Protocol:

1. Inoculate 5 mL LB media with a single colony from plate.
2. Grow cells overnight at 30 °C, 225 rpm.

## Day 3: Probe cells and extract RNA (~10 hours)

Materials:

- LB media (50 mL)
- PBS buffer (100 mL)

- 50 mM ligand stock in DMSO (20 uL per sample)
- Sterile 6-well plastic plate (1 well per sample)
- lysozyme (60 mg)
- 1M Tris pH 7 (180 uL)
- 0.5 M EDTA (120 uL)
- Chloroform
- 100% ethanol
- Qiagen RNAeasy Plus Universal Mini kit

#### Protocol:

3. Inoculate 60 mL LB media with 3 mL cultured cells and grow cells to 0.50 OD<sub>600</sub>.

**NOTE: Expect cell growth to take ~2.5 hours.**

Time point	#1 OD <sub>600</sub>	#2 OD <sub>600</sub>	#3 OD <sub>600</sub>
	-	-	-

4. Buffer exchange cells.

- Gently pellet cells at 3200g for 10 min at 4 °C.
- Carefully remove supernatant.
- Wash cell pellet with 50 mL PBS buffer, repeat centrifugation, and carefully remove supernatant.
- Resuspend cells in 50 mL PBS buffer.

5. Probe cells with functionalized ligands.

- Combine 4,980 uL of cultured media with 20 uL of photo-reactive ligand (50 mM stock in DMSO) and pipette mix. Avoid prolonged light exposure.

**NOTE: final [ligand] is 200 uM. 250X dilution.**

- Incubate ligand-treated cells for 30 minutes at 30 °C and 225 rpm.

**NOTE: Keep solution in dark.**

6. Cross-link functionalized ligands to RNA.

- Transfer 5 mL cell solution to a sterile 6-well plastic plate.
- Place uncover plate and place on ice
- Irradiate with 365 nm light source and cross-link with 3 J/cm<sup>2</sup>.

**NOTE: This is nine minutes (or three max energy hits) 5-6 cm from the light source in a UVP crosslinker.**

7. Create cell pellet.

- Centrifuge to pellet the cells. 8,000xg for 10 min. at 4 °C.

- Discard supernatant.

8. Prepare fresh lysis buffer.

- (30 mM Tris pH 7, 10 mM EDTA, 10 mg/mL lysozyme):

Reagent	Amount
lysozyme	50 mg
1M Tris pH 7	150 uL
0.5M EDTA	100 uL
RNase inhibitor murine	50 uL
NF H2O	4700 uL
<b>TOTAL</b>	<b>5,000 uL</b>

9. Lysis cell pellet with lysozyme buffer.

- Add 200 uL of buffer lysis to each cell pellet.
- Incubate for 30 min. at RT.

10. Extract and purify total RNA from cells.

- Use 1000 uL of Tri-reagent
- Mix aggressively to aid in lysis
- Lyse for 15 minutes

**OPTIONAL PAUSE: store at 4 °C**

- Add 250 uL of chloroform
- Vortex and incubate at RT for 3-5 minutes
- Centrifuge at 15000xg for 15 minutes at 4 °C

**NOTE: Keep sample on ice until aqueous layer is recovered**

- Add 100% ethanol to each recovered aqueous layer in equal parts (1:1)
- Follow [Monarch total RNA miniprep](#) kit instructions. Start at Part 2 step 4.
- Elute from column with 51.5 uL of NFW

**NOTE: Yields of ~50-60 ug total RNA have been achieved per 5 mL culture probed.**

11. Check RNA purity.df

- Measure absorbance ratio (A260/280 and A260/230) on nanodrop.

Sample ID	Replicate	Ligand	Yield (µg)	[RNA] (ng/µL)	A260/280	A260/230	Vol NFW to 500 ng/µL (µL)

NOTES:

Storage:

- RNA can be stored overnight at 4 °C for ~ 1 week.
- RNA can be flash-frozen in aliquots and stored at -80 °C for ~1 year.

Next Steps:

- Partial RNA Fragmentation